FINAL REPORT

Sino-implant (II) Initiative: 2012 Quality Assurance Evaluation

December 15, 2012

Submitted by:

FHI 360

Markus Steiner, Ph.D. Sino-implant (II) Project Director P.O. Box 13950 Research Triangle Park, NC 27709

Tel: +1.919.544.7040 FAX: +1.919.544.7261

Final Report

Sino-implant (II): 2012 Quality Assurance Evaluation

FHI 360

David W. Jenkins, Ph.D. Associate Scientist II

Product Quality and Compliance Department 2810 Meridian Parkway, Suite 133 Durham, NC 27713

E-mail: djenkins@fhi360.org Phone: 919-544-7040 ext. 11617

Fax: 919-544-5849

Derek H. Owen, Ph.D. Scientist I

Clinical Sciences 2224 E. NC Hwy. 54, Durham NC 27713

E-mail: dowen@fhi360.org

Telephone: 919-544-7040 ext. 11168

Fax: 919-544-7261

Markus Steiner, Ph.D. Sino-implant (II) Project Director

Clinical Sciences

2224 E. NC Hwy. 54, Durham NC 27713

E-mail: msteiner@fhi360.org

Telephone: 919-544-7040 ext. 11346

Fax: 919-544-7261

Table of Contents

Final Report	i
Executive Summary	iv
Introduction	1
Section I: Quality Assurance Monitoring of Sino-implant (II)	2
Brief Description of Sino-implant (II) Commercial Lot Release Testing	
Results of Sino-implant (II) Commercial Lot Release Verification	
Conclusion of 2012 Quality Assurance Monitoring of Commercial Lots of Sino-implant	
Section II: Annual Quality Assurance Evaluation of Sino-implant (II)	4
Levonorgestrel: Active Pharmaceutical Ingredient (API) Quality Assurance Evaluation	
Levonorgestrel CP Evaluation	4
Levonorgestrel BP Evaluation	5
Conclusion of Levonorgestrel: Active Ingredient Quality Assurance Evaluation	6
Sino-implant (II): Final Product Quality Assurance Evaluation	7
Sino-implant (II) Sample Lot Release Verification	
Ethylene Oxide Residuals Evaluation	8
Inorganic Impurities Evaluation	9
Residual Solvents Evaluation	13
Bacterial Endotoxin Evaluation	14
Cytotoxicity Evaluation	15
Sino-implant (II): Packaging Material Evaluation	15
Packaging Physicochemical Evaluation	16
Package Integrity Evaluation	16
Section III: Summary of Results of the 2012 Comprehensive Quality Assurance Evaluation	18

Tables

Table 1. Summary of the quality assurance evaluation of Sino-implant (II)	V
Table 2. Sino-implant (II) commercial lot release verification test results	3
Table 3. Levonorgestrel: lot release verification test results analyzed using the Chinese	
Pharmacopeia monograph	5
Table 4. Levonorgestrel: lot release verification test results analyzed using the BP monograph.	5
Table 5. Sino-implant (II) sample lot release verification results	7
Table 6. Ethylene oxide residuals test results of Sino-implant (II)	8
Table 7. ICP-MS test results of Sino-implant (II): USP inorganic impurities (Heavy Metals)	11
Table 8. ICP-MS test results of Sino-implant (II): additional inorganic impurities	12
Table 9. Residual solvents test results of Sino-implant (II)	4
Table 10. Limulus Amebocyte Lysate Test (Kinetic Turbidimetric Technique) test results of	
Sino-implant (II).	4
Table 11. MEM elution test results of Sino-implant (II).	15
Table 12. Physicochemical test results of Sino-implant (II) packaging film	6
Table 13. Bubble emission test results of Sino-implant (II).	17
Table 14. Comprehensive tabulated summary of the quality assurance results of Sino-implant	
(II) tests conducted	8

iii

Executive Summary

A quality assurance evaluation and monitoring program was implemented by FHI 360 to verify that Sino-implant (II) meets lot release specifications for the product. This report expands on the results obtained after the quality assurance evaluation and monitoring program conducted from 2008 to 2011.

Commercial lot release testing verification was conducted for all eight lots of Sino-implant (II) that were shipped to countries as part of the Sino-implant (II) introduction in developing countries. In all instances results were in compliance with Sino-implant (II) lot release specifications (Table 1).

A more extensive battery of tests was conducted to evaluate if Sino-implant (II) meets international quality standards. The battery of tests selected was based on standards from the United States Pharmacopeia (USP), British Pharmacopeia (BP), International Organization for Standardization (ISO) and the American Society for Testing and Materials (ASTM). Sino-implant (II) samples tested met the quality requirements for all tests conducted (Table 1).

We conclude that Shanghai Dahua Pharmaceutical Co., Ltd continues to demonstrate the ability to consistently produce a contraceptive implant that meets international quality standards.

Table 1. Summary of the quality assurance evaluation of Sino-implant (II)

Test		Lot #	Results			
Quality Assurance Monitoring of Sino-implant (II)						
Sino-implant (II): Commercial Lot Release Testing Verification ^a		07112001 05162012 06082012 06272012 07132012 08032012 09152012 10072012	Met requirements			
Annual Quality Assurance Evaluation of Sino-implant (II)						
Levonorgestrel: Active Ingredient Quality Assurance Evaluation	Levonorgestrel: China Pharmacopeia Evaluation ^a	ZQ20120002 ZQ20110016 ZQ20110022	Met requirements			
	Levonorgestrel: British Pharmacopeia Evaluation ^e	2020110022	Met requirements			
Sino-implant(II): Final Product Quality Assurance Evaluation	Sample Lot Release Verification ^a	09022012 24022012 10032012	Met requirements			
	Ethylene Oxide Residuals Evaluation ^b	10032012	Met requirements			
	Inorganic Impurities Evaluation ^c		Met requirements			
	Residual Solvents Evaluation ^d		Met requirements			
	Bacterial Endotoxin Evaluation ^b		Met requirements			
	Cytotoxicity Evaluation ^b		Met requirements			
Sino-implant(II): Packaging Material Evaluation	Packaging Physicochemical Evaluation ^b	11-071PA1-11-OP-22-1 11-071P-11-8-23-1 110117	Met requirements			
	Sino-implant (II): Package Integrity Evaluation ^b	09022012 24022012 10032012	Met requirements			

^a Conducted by SGS Life Sciences Division, Shanghai, China ^b Conducted by Nelson Laboratories, Inc., Salt Lake City, UT ^c Conducted by EAG-Life Sciences, Maryland Heights, MO ^d Conducted by Irvine Pharmaceutical Sciences, Irvine, CA ^e Conducted by Lancaster Laboratories, Lancaster, PA